



Aequitas-Medical

Kerkstraat 68 TEL +32 2 378.29.32
1601 Ruisbroek FAX +32 2 331.56.73
BELGIUM
www.aequitas-medical.be
info@aequitas-medical.be

 **FloShieldAir**[™]
Superior Laparoscopic Vision

FloShield Air[™]

Active Lens Protection

A cost effective, all-inclusive solution for maintaining laparoscopic clarity

FloShield Air delivers continuous defogging during laparoscopic surgery while dramatically reducing the need to remove the scope for cleaning.

Continuously Defogs


Pushes Away Smoke & Debris

Competitively Priced

Dramatically Reduces Scope Removals



The FloShield Difference



Active Lens Protection



Defogs the Lens In Situ

Vortex Barrier Technology provides a continuous flow of dry CO₂ gas over the tip of the scope, keeping the lens consistently clear of condensation and fog... *Without removing the scope from the surgical site.*

Creates a CO₂ Barrier

The FloShield Air sheath directs CO₂ from the surgical insufflator to the tip of the laparoscope, creating a barrier of dry gas over the lens that shields it from condensation, smoke and debris.

FloShield Air actively defogs while creating a CO₂ barrier that shields the lens from smoke and debris.



FloShield Air Fits

Most 5mm or 10mm scopes
0°, 30°, and 45° angled lens

FloShield Air Provides Enhanced Clinical Benefits that Deliver Better Value than Clarify™

	FloShield Air	Clarify™ (formerly D-HELP®)
Instantly Defogs	★★★★	★
Continuously Defogs	★★★★	★
Protective CO ₂ Barrier	yes	no
Dramatically Reduces Scope Removals	yes	no
Eliminates Need to Heat Lens	yes	no
Battery Disposal	n/a	yes
Competitively Priced	★★★★	★★★★

Clarify (D-HELP) is a registered trademark of Covidien.

www.floshield.com