







## Cervix Adaptor

The Cervix Adaptor for pertubation, hystero-salpingo-graphy, and hydro-pertubation. The cervical canal can be sealed easily and without any pain. The plastic adaptor, 3 sizes available, adapts itself to any cervical configuration without altering the topography of the interior genitals

	Art. No.:	
	<b>1235 A</b> <b>1235 B</b> <b>1235 C</b>	Cervix Adaptor, disposable, for insufflation and HSG, sterile individual wrap, tube length: 330 mm, <b>10 pcs p.pck.</b> Cup diameter 25 mm Cup diameter 27 mm Cup diameter 30 mm
	<b>1240 A</b> <b>1240 B</b> <b>1240 C</b>	Cervix Adaptor, disposable, for hydro-pertubation, sterile individual wrap, tube length: 330 mm, <b>10 pcs p.pck.</b> cup diameter 25 mm cup diameter 27 mm cup diameter 30 mm
	<b>1245 A</b> <b>1245 B</b> <b>1245 C</b>	Cervix Adaptor, disposable, for insemination and drug delivery at the cervix, sterile individual wrap, tube length: 220 mm, <b>10 pcs p.pck.</b> cup diameter 25 mm cup diameter 27 mm cup diameter 30 mm
	<b>1246 A</b> <b>1246 B</b> <b>1246 C</b>	Cervix Cups, disposable, for artificial insemination and drug delivery at the cervix, sterile individual wrap, threath length: 330 mm, <b>10 pcs p.pck.</b> cup diameter 25 mm cup diameter 27 mm cup diameter 30 mm
	<b>1260</b>	Grasping Forceps for Cervix Adaptor, made of stainless steel
	<b>1261</b>	Vacuum Handpump, with control manometer and accessories: silicone connection tubing, sterilizable secretion trap

## Uterus Manipulator



- Ergonomic handle for left and right hand.
- One hand operation.
- Stainless steel version, incl. adaptors.
- Without screws und nuts.
- Automatically ratchet.

**Uterus Manipulator,  
inclusive:**

1267 S

**Working length Diameter**

**Cone, with short thorn**

1267-1



60 mm

25 mm

**Cone, with long thorn**

1267-2



70 mm

25 mm

**Uterus Obturator, thick, extra long**

1267-6




100mm

10 mm

**Hook Forceps**

1265



*Aequitas-Medical*

Distributor Belgium :

Aequitas-Medical  
Kerkstraat 68 – 1601 Ruisbroek,  
Belgium  
TEL +32 2 378 29 32  
FAX +32 2 331 56 73

[info@aequitas-medical.be](mailto:info@aequitas-medical.be) / [www.aequitas-medical.be](http://www.aequitas-medical.be)