

## Dispo seal Morcellator

Continuous Extraction Morcellator

15 mm für / for

MORCDRIVE™ / POWERDRIVE

- Multiple Use Morcellator with disposable seal.
- Sterile packed disposable seal for an easy replacement of the most stressed part of the morcellation system.
- The loss of CO<sub>2</sub> is prevented by a silicone valve in the disposeal, which allows for the use of the grasping instrument.
- The morcellator allows the surgeon to directly extract morcellated tissue without the need to remove the complete device. This simple to use device will increase the patient safety and decrease the operating time of the procedure.
- When morcellating soft tissues the structure can vary quite considerably from patient to patient.
- Whether the tissue mass is soft and compliant or fibrosus and resilient, our cutting tubes have been specially designed to make easy work of the procedure.
- Uncomplicated and time-saving morcellation.
- Riskless intra-abdominal, intrathoracal morcellation, even of larger tissues.
- Sharp and resistant edge of tube.
- Cutting tube with great stability, made of high-quality stainless steel.
- Working length 240 mm.



Morc Drive / Power Drive with inserted Dispo seal Morcellator



*Aequitas-Medical*

Distributor Belgium : Aequitas-Medical  
Kerkstraat 68 – 1601 Ruisbroek,  
Belgium  
TEL +32 2 378 29 32  
FAX +32 2 331 56 73

[info@aequitas-medical.be](mailto:info@aequitas-medical.be) / [www.aequitas-medical.be](http://www.aequitas-medical.be)



The dispo seal is simply screwed into the reusable morcellator.

Diameter	Art. No.	Morc Drive / Power Drive
----------	----------	--------------------------



15 mm	7689 W1D
15 mm	7689 M1D
	7689 DKE5

Dispo seal Morcellator  
Continuous Extraction Morcellator  
for 10 mm Instruments,  
incl. disposable seal 7689 DKE  
(5 pcs.p.pck.)

Wave Edge



Serrated Edge

Disposable seal, sterile coated, 5 pcs.p.pck

### Power Drive™ for Dispo seal Morcellator



**Power Drive™ set**

with  
Gearbox  
f. Macro Morc

**7688 PD1**

7688 PDU  
7688 PDC  
7688 PDMO  
7688 PDG1

**Power-Drive**  
consisting of:

Controlling unit  
Connecting cable, 3 m  
Drive Motor  
Gearbox f. Macro Morc  
Mains cable 2,5 m