

**Twin Cut™**

Continuous Extraction Morcellator  
15 mm for  
**MORCDRIVE™ / POWERDRIVE**

- Counter Rotating Twin-Cut Blades to reduce tissue rotation.
- The morcellator allows the surgeon to directly extract morcellated tissue without the need to remove the complete device. This simple to use device will increase the patient safety and decrease the operating time of the procedure.
- The Continuous extraction morcellator blades come as a single use product. When morcellating soft tissues the structure can vary quite considerably from patient to patient.
- Whether the tissue mass is soft and compliant or fibrous and resilient, our cutting tubes have been specially designed to make easy work of the procedure.
- Uncomplicated and time-saving morcellation.
- Riskless intra-abdominal, intra-thoracal morcellation, even of larger tissues.
- Sharp and resistant edge of tube.
- Cutting tube with great stability, made of high-quality stainless steel.
- Working length 240 mm.



Morc Drive mini / Power Drive with inserted Twin-Cut Morcellator with counter rotating cutting edges



*Aequitas-Medical*

Distributor Belgium :

Aequitas-Medical

Kerkstraat 68 – 1601 Ruisbroek,

Belgium

TEL +32 2 378 29 32

FAX +32 2 331 56 73

[info@aequitas-medical.be](mailto:info@aequitas-medical.be) / [www.aequitas-medical.be](http://www.aequitas-medical.be)

Diameter

Art. No

Morc Drive mini / Power Drive



Twin Cut, Morcellators with counter rotating cutting edges. Continuous Extraction Morcellator for 10 mm Instruments, 5 pcs.p.pck.



15 mm

7689 W1DWD5

Wave Edge



15 mm

7689 M1DSD5

Serrated Edge



**Morc Drive™ mini set**

with Gearbox f. Twin Cut



**Power Drive™ set**

with Gearbox f. Twin Cut



**7688 MINT**

- 7688 MINM
- 7688 PDGD
- 7688 MINLW

**MORCDRIVE™ mini set**

consisting of:  
 MORCDRIVE™ , 2 pcs.  
 Gearbox double f. Twin Cut  
 Mini-Charger with medical external power supply 220 - 240 V

**7688 PD1T**

- 7688 PDU
- 7688 PDC
- 7688 PDMO
- 7688 PDGD

**Power-Drive:**

consisting of  
 Controlling unit  
 Connecting cable, 3 m  
 Drive Motor  
 Gearbox double f. Twin Cut  
 Mains cable 2,5 m