




SUS Jet – S064 Suction Irrigator



- Fast and effective irrigation
- 50% less fluid use – reduced suction time
- Clean wash (No cross contamination)
- Smooth and leak-free valves
- Simultaneous suction and irrigation

Jet flow or...



... High flow

Laparoscopy - Suction/Irrigator

Part financed by
EU European
Regional
Development Fund




Aequitas-Medical

Aequitas-Medical
Kerkstraat 68
1601 Ruisbroek
BELGIUM

TEL +32 2 378.29.32
FAX +32 2 331.56.73
www.aequitas-medical.be
info@aequitas-medical.be


Features & Benefits



SUS Suction Irrigator	Trumpet Valve Suction Irrigator
<ul style="list-style-type: none"> Pistol-grip & dual tubing: comfortable wrist/hand position, better control 	<ul style="list-style-type: none"> Bulky, heavy handles with Inline-grip: uncomfortable & poor control
<ul style="list-style-type: none"> Choice of 2 flows - Jet & High Flow Faster, more effective irrigation Better targeting 50% less fluid use Ability to simultaneously aspirate and irrigate 	<ul style="list-style-type: none"> Single Channel for suction and irrigation contaminating irrigation fluid and providing poor visibility Large flow / low pressure
<ul style="list-style-type: none"> Dual Channels providing clean wash Irrigation and waste fluids are separated allowing clearer vision 	<ul style="list-style-type: none"> Irrigation fluid contamination by single tube
<ul style="list-style-type: none"> Soft seal valves: smooth, no leaks 	<ul style="list-style-type: none"> Valves leak / stick. Uncomfortable to operate
<ul style="list-style-type: none"> Clean room manufactured, supplied sterile 	<ul style="list-style-type: none"> Hard to clean, valves need stripping before washing. Cleanliness issue remains

S063 - Irrigation Pump

The suction irrigator package is completed with a reusable irrigation pump



- Quick and easy setup
- Constant flow
- Adjustable pressure
- Clamp to drip stand
- No need for cleaning
- No pressure cuff required

Laparoscopy - Suction/Irrigator

Part financed by
EU European
Regional
Development Fund



Aequitas-Medical

Aequitas-Medical
Kerkstraat 68
1601 Ruisbroek
BELGIUM

TEL +32 2 378.29.32
FAX +32 2 331.56.73
www.aequitas-medical.be
info@aequitas-medical.be