

Single Cut™

Continuous Extraction Morcellator

10,15 mm for

MORCDRIVE™ mini/ POWERDRIVE

- The morcellator allows the surgeon to directly extract morcellated tissue without the need to remove the complete device. This simple to use device will increase the patient safety and decrease the operating time of the procedure.
- The Continuous extraction morcellator blades come as a single use product. When morcellating soft tissues the structure can vary quite considerably from patient to patient.
- Whether the tissue mass is soft and compliant or fibrous and resilient, our cutting tubes have been specially designed to make easy work of the procedure.
- Uncomplicated and time-saving morcellation.
- Riskless intra-abdominal, intra-thoracal morcellation, even of larger tissues.
- Sharp and resistant edge of tube.
- Cutting tube with great stability, made of high-quality stainless steel.
- Working length 240 mm.



Morc Drive mini / Power Drive with inserted Single Cut Morcellator

*Aequitas-Medical*

Distributor Belgium :

Aequitas-Medical

Kerkstraat 68 – 1601 Ruisbroek,

Belgium

TEL +32 2 378 29 32

FAX +32 2 331 56 73

info@aequitas-medical.be / www.aequitas-medical.be

Diameter

Art. No

Morc Drive mini / Power Drive



Single Cut, Continuous Extraction Morcellator,
5 pcs.p.pck.



10 mm
15 mm

7689 W6DWE
7689 W1DWE

Wave Edge
Wave Edge



10 mm
15 mm

7689 M6DSE
7689 M1DSE

Serrated Edge
Serrated Edge

**Morc
Drive™
mini set**

with
Gearbox f.
Single Cut



**Power
Drive™
set**

with
Gearbox f.
Single Cut



7688 MINTS

7688 MINM

7688 PDG

7688 MINLW

MORCDRIVE™ mini set

consisting of:

MORCDRIVE™, 2 pcs.

Gearbox f. Single Cut

Mini-Charger with medical external
power supply 220 - 240 V

7688 PD1S

7688 PDU

7688 PDC

7688 PDMO

7688 PDG

Power-Drive:

consisting of

Controlling unit

Connecting cable, 3 m

Drive Motor

Gearbox f. Single Cut

Mains cable 2,5 m